


ALIX Rabbit mAb

ABZ1974

Product Overview

Name	ALIX Rabbit mAb
Accession(Primary)	Q8WUM4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 96 kDa; Observed MW: 96 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q8WUM4)

Gene Name PDCD6IP**About Gene Symbol**

Western blot analysis of ALIX in CHO-K1, rat Brain, Jurkat lysates using ALIX antibody.

