

## CD276 Rabbit mAb


ABZ1978

### Product Overview

|                             |                                                                          |
|-----------------------------|--------------------------------------------------------------------------|
| <b>Name</b>                 | CD276 Rabbit mAb                                                         |
| <b>Accession(Primary)</b>   | Q8VE98                                                                   |
| <b>Host Species</b>         | Rabbit                                                                   |
| <b>Clonality</b>            | Monoclonal                                                               |
| <b>Isotype</b>              | IgG                                                                      |
| <b>Species Reactivity</b>   | Human, Mouse                                                             |
| <b>Calculated MW</b>        | Calculated MW: 34 kDa; Observed MW: 90-110 kDa                           |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; ICC/IF-1:500-1:1000;<br>FC-1:1000-1:5000; IP-1:20-1:100 |

### Target information (Q8VE98)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | Cd276 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of CD276 in HepG2 lysates  
using CD276 antibody.