


LRRC15 Rabbit mAb

ABZ1979

Product Overview

Name	LRRC15 Rabbit mAb
Accession(Primary)	Q8TF66
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64 kDa; Observed MW: 64 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:100

Target information (Q8TF66)

Gene Name LRRC15**About Gene Symbol**

Western blot analysis of LRRC15 in rat Brain, C6 lysates using LRRC15 antibody.

