


KMT4 Rabbit mAb

ABZ1980


Product Overview

Name	KMT4 Rabbit mAb
Accession(Primary)	Q8TEK3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 165 kDa; Observed MW: 185 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:100-1:200

Target information (Q8TEK3)

Gene Name DOT1L**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded mouse testis using KMT4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of KMT4/Dot1L in HeLa, A549, HL-60 lysates using KMT4/Dot1L antibody.

Copyright© 2024 zen-bioscience Limited.