


UBA3 Rabbit mAb

ABZ1983


Product Overview

Name	UBA3 Rabbit mAb
Accession(Primary)	Q8TBC4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 52 kDa
Recommended Dilution	WB-1:2000-1:5000

Target information (Q8TBC4)

Gene Name UBA3**About Gene Symbol**

Western blot analysis of UBA3 in Raw264.7 lysates using NEDD8 Activating Enzyme E1C antibody.



Western blot analysis of UBA3 in HeLa, A549, HL-60 lysates using UBA3 antibody.