


## SEH1L Rabbit mAb

ABZ1984

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | SEH1L Rabbit mAb                           |
| <b>Accession(Primary)</b>   | Q8R2U0                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse                               |
| <b>Calculated MW</b>        | Calculated MW: 40 kDa; Observed MW: 40 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IP-1:20-1:50              |

### Target information (Q8R2U0)

**Gene Name** Seh1L**About Gene Symbol**

Western blot analysis of SEH1L in K562, C6, 3T3 lysates using SEH1L antibody.

