

HOPX Rabbit mAb

ABZ1985


Product Overview

Name	HOPX Rabbit mAb
Accession(Primary)	Q8R1H0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 8 kDa; Observed MW: 12 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (Q8R1H0)

Gene Name Hopx

About Gene Symbol



Western blot analysis of HOPX/HOD in HeLa lysates using HOPX antibody.

