

Nesprin1 Rabbit mAb


ABZ1986

Product Overview

Name	Nesprin1 Rabbit mAb
Accession(Primary)	Q8NF91
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 1011 kDa; Observed MW: 112 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200

Target information (Q8NF91)

Gene Name	SYNE1
About Gene Symbol	



Western blot analysis of Nesprin 1 in MCF-7, C6, 3T3, HeLa lysates using Nesprin1 antibody.