


Transmembrane Protein 161B Rabbit mAb

ABZ1987

Product Overview

Name	Transmembrane Protein 161B Rabbit mAb
Accession(Primary)	Q8NDZ6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 55 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000; IP-1:20-1:50

Target information (Q8NDZ6)

Gene Name TMEM161B**About Gene Symbol**

Western blot analysis of TMEM161B in K562, HeLa lysates using Transmembrane Protein 161B antibody.

