

HJURP Rabbit mAb

ABZ1988


Product Overview

Name	HJURP Rabbit mAb
Accession(Primary)	Q8NCD3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100


Target information (Q8NCD3)

Gene Name HJURP

About Gene Symbol



Western blot analysis of HJURP in HeLa, CHO-K1 lysates using HJURP antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HJURP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.