

FANCB Rabbit mAb


ABZ1990

Product Overview

Name	FANCB Rabbit mAb
Accession(Primary)	Q8NB91
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 98 kDa; Observed MW: 98 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; FC-1:50-1:100; IP-1:20-1:50

Target information (Q8NB91)

Gene Name	FANCB
About Gene Symbol	



Western blot analysis of FANCB in 3T3, HeLa
lysates using FANCB antibody.