


FNDC5 Rabbit mAb

ABZ1991

Product Overview

Name	FNDC5 Rabbit mAb
Accession(Primary)	Q8NAU1
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q8NAU1)

Gene Name FNDC5**About Gene Symbol**

Western blot analysis of FNDC5 in K562, C6, 3T3, HeLa lysates using FNDC5 antibody.

