

SCAI Rabbit mAb


ABZ1992

Product Overview


Name	SCAI Rabbit mAb
Accession(Primary)	Q8N9R8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 60 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:200-1:500; FC-1:100-1:200

Target information (Q8N9R8)

Gene Name	SCAI
About Gene Symbol	



Western blot analysis of SCAI in C6, 3T3 lysates using SCAI antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using SCAI antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.