

SIRT6 Rabbit mAb


ABZ1993

Product Overview


Name	SIRT6 Rabbit mAb
Accession(Primary)	Q8N6T7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 39 kDa; Observed MW: 39 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:500-1:1000; IP-1:20-1:50

Target information (Q8N6T7)

Gene Name	SIRT6
About Gene Symbol	



Western blot analysis of SIRT6 in 3T3, C6 lysates using SIRT6 antibody.



Western blot analysis of SIRT6 in 293T, Jurkat lysates using SIRT6 antibody.