

DBC 1 Rabbit mAb


ABZ1994

Product Overview

Name	DBC 1 Rabbit mAb
Accession(Primary)	Q8N163
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 103 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:200-1:500; FC-1:200-1:1000; IP-1:20-1:50

Target information (Q8N163)

Gene Name	CCAR2
About Gene Symbol	



Western blot analysis of DBC1 in 293, Jurkat
lysates using DBC 1 antibody.