

CD226 Rabbit mAb


ABZ1995

Product Overview

Name	CD226 Rabbit mAb
Accession(Primary)	Q8K4F0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 38 kDa; Observed MW: 55-65 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (Q8K4F0)

Gene Name	Cd226
About Gene Symbol	



Western blot analysis of CD226 in mouse spleen
lysates using CD226 antibody.