

## PELP1 Rabbit mAb


ABZ1996

### Product Overview


|                             |                                                                           |
|-----------------------------|---------------------------------------------------------------------------|
| <b>Name</b>                 | PELP1 Rabbit mAb                                                          |
| <b>Accession(Primary)</b>   | Q8IZL8                                                                    |
| <b>Host Species</b>         | Rabbit                                                                    |
| <b>Clonality</b>            | Monoclonal                                                                |
| <b>Isotype</b>              | IgG                                                                       |
| <b>Species Reactivity</b>   | Human                                                                     |
| <b>Calculated MW</b>        | Calculated MW: 120 kDa; Observed MW: 160 kDa                              |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:100-1:200;<br>ICC/IF-1:100-1:200; FC-1:200-1:500 |

### Target information (Q8IZL8)


|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | PELP1                                                                                                                                                                                                                                                           |
| <b>About Gene Symbol</b> | PELP1 (Proline, Glutamate And Leucine Rich Protein 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mec |



Western blot analysis of PELP1 in HeLa, A549, MCF-7, 293T lysates using PELP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PELP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of PELP1 (green)  
in hela using PELP1 antibody, and DAPI(blue)