

# RAB45 Rabbit mAb

ABZ1997


## Product Overview

|                             |                                                  |
|-----------------------------|--------------------------------------------------|
| <b>Name</b>                 | RAB45 Rabbit mAb                                 |
| <b>Accession(Primary)</b>   | Q8IZ41                                           |
| <b>Host Species</b>         | Rabbit                                           |
| <b>Clonality</b>            | Monoclonal                                       |
| <b>Isotype</b>              | IgG                                              |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                |
| <b>Calculated MW</b>        | Calculated MW: 83 kDa; Observed MW: 83 kDa       |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:20-1:50 |


## Target information (Q8IZ41)

**Gene Name** RASEF


**About Gene Symbol**



Western blot analysis of RAB45 in mouse heart lysates using Rab45 antibody.



Western blot analysis of RAB45 in Jurkat lysates using RAB45 antibody.



Immunocytochemistry analysis of RAB45 (green) in HeLa using RAB45 antibody, and DAPI (blue)