


MSRB3 Rabbit mAb

ABZ1998


Product Overview

Name	MSRB3 Rabbit mAb
Accession(Primary)	Q8IXL7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q8IXL7)

Gene Name MSRB3**About Gene Symbol**

Western blot analysis of MSRB3 in 3T3 lysates using MSRB3 antibody.



Western blot analysis of MSRB3 in C6 lysates
using MSRB3 antibody.