

SIRT2 Rabbit mAb

ABZ1999


Product Overview

Name	SIRT2 Rabbit mAb
Accession(Primary)	Q8IXJ6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q8IXJ6)

Gene Name SIRT2

About Gene Symbol



Western blot analysis of SIRT2 in K562, HeLa lysates using SIRT2 antibody.

