


## SIRT2 Rabbit mAb

ABZ2000


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | SIRT2 Rabbit mAb                           |
| <b>Accession(Primary)</b>   | Q8IXJ6                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 43 kDa; Observed MW: 39 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-1:50-1:100; IP-1:20   |


### Target information (Q8IXJ6)

**Gene Name** SIRT2**About Gene Symbol**

Western blot analysis of SIRT2 in mouse brain,  
U87-MG lysates using SIRT2 antibody.



Western blot analysis of SIRT2 in rat Brain lysates using SIRT2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using SIRT2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.