


TIAM2 Rabbit mAb

ABZ2001

Product Overview

Name	TIAM2 Rabbit mAb
Accession(Primary)	Q8IVF5
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 190 kDa; Observed MW: 190 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200

Target information (Q8IVF5)

Gene Name TIAM2**About Gene Symbol**

Western blot analysis of TIAM2 in K562, 293 lysates using TIAM2 antibody.

