


TRIF Rabbit mAb


ABZ2002

Product Overview

Name	TRIF Rabbit mAb
Accession(Primary)	Q8IUC6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 76 kDa; Observed MW: 98 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50

Target information (Q8IUC6)

Gene Name TICAM1**About Gene Symbol**



Western blot analysis of TRIF in K562, 3T3, HeLa lysates using TRIF antibody.