

Caspase 9 Rabbit mAb

ABZ2003


Product Overview

Name	Caspase 9 Rabbit mAb
Accession(Primary)	Q8C3Q9
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20-1:50


Target information (Q8C3Q9)

Gene Name Casp9

About Gene Symbol



Western blot analysis of Caspase9 in 3T3 lysates using Caspase 9 antibody.



Western blot analysis of Caspase9 in C6 lysates
using Caspase9 antibody.