

LYVE1 Rabbit mAb


ABZ2004

Product Overview


Name	LYVE1 Rabbit mAb
Accession(Primary)	Q8BHC0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35-70 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:1000-1:5000; FC-1:50-1:100; IP-1:20-1:50

Target information (Q8BHC0)

Gene Name	Lyve1
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse lung using LYVE1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of LYVE1 in mouse Lung lysates using LYVE1 antibody.