

Homer1 Rabbit mAb


ABZ2005

Product Overview


Name	Homer1 Rabbit mAb
Accession(Primary)	Q86YM7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500; IP-1:50-1:100

Target information (Q86YM7)

Gene Name	HOMER1
About Gene Symbol	HOMER1 (Homer Scaffold Protein 1) is a signaling scaffold protein involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Homer1 in mouse brain lysates using Homer1 antibody.



Western blot analysis of Homer1 in rat Brain lysates using Homer1 antibody.