

NIT1 Rabbit mAb

ABZ2006

Product Overview

Name	NIT1 Rabbit mAb
Accession(Primary)	Q86X76
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:20

Target information (Q86X76)

Gene Name NIT1**About Gene Symbol**

