

STING Rabbit mAb


ABZ2007

Product Overview


Name	STING Rabbit mAb
Accession(Primary)	Q86WV6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 42 kDa; Observed MW: 42 kDa
Recommended Dilution	WB-1:1000-1:2000; IHC-1:100-1:1000; ICC/IF-1:500-1:2000; FC-1:20; IP-1:20-1:50

Target information (Q86WV6)

Gene Name	STING1
About Gene Symbol	STING1 (Stimulator Of Interferon Response CGAMP Interactor 1) is an innate immune sensor involved in JAK-STAT/cytokine signaling and innate immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely



Western blot analysis of TMEM173 in K562, rat Brain lysates using Transmembrane Protein 173 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using TMEM173 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.