

Aly Rabbit mAb


ABZ2009

Product Overview


Name	Aly Rabbit mAb
Accession(Primary)	Q86V81
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 27 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:500; ICC/IF-1:100-1:200; FC-1:50-1:100; IP-1:20-1:50

Target information (Q86V81)

Gene Name	ALYREF
About Gene Symbol	



Western blot analysis of Aly/Ref in rat brain, rat heart lysates using Aly antibody.



Western blot analysis of Aly/Ref in HepG2, mouse liver lysates using Aly/Ref antibody.