

NgR3 Rabbit mAb


ABZ2010

Product Overview


Name	NgR3 Rabbit mAb
Accession(Primary)	Q86UN2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 49 kDa; Observed MW: 49 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:500; ICC/IF-1:50-1:100

Target information (Q86UN2)


Gene Name	RTN4RL1
About Gene Symbol	



Western blot analysis of NgR3 in rat brain, mouse brain lysates using NgR3 antibody.



Western blot analysis of NgR3 in A549, HL-60, U2OS lysates using NgR3 antibody.



Immunocytochemistry analysis of NgR3 (green) in Hela using NgR3 antibody, and DAPI (blue)