

LYRIC/AEG1 Rabbit mAb


ABZ2011

Product Overview


Name	LYRIC/AEG1 Rabbit mAb
Accession(Primary)	Q86UE4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 64 kDa; Observed MW: 75 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:2000; ICC/IF-1:100; FC-1:200-1:500; IP-1:20-1:50

Target information (Q86UE4)


Gene Name	MTDH
About Gene Symbol	



Immunocytochemistry analysis of LYRIC (green) in Hela using LYRIC antibody, and DAPI (blue).



Western blot analysis of LYRIC/AEG1 in Hela, MCF-7, HepG2, U251, HL-60 lysates using LYRIC/AEG1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using LYRIC/AEG1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.