

## LYRIC/AEG1 Rabbit mAb


ABZ2012

### Product Overview


|                             |                                                            |
|-----------------------------|------------------------------------------------------------|
| <b>Name</b>                 | LYRIC/AEG1 Rabbit mAb                                      |
| <b>Accession(Primary)</b>   | Q86UE4                                                     |
| <b>Host Species</b>         | Rabbit                                                     |
| <b>Clonality</b>            | Monoclonal                                                 |
| <b>Isotype</b>              | IgG                                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                          |
| <b>Calculated MW</b>        | Calculated MW: 64 kDa; Observed MW: 75 kDa                 |
| <b>Recommended Dilution</b> | WB-1:2000-1:5000; IHC-P-1:100-1:200;<br>ICC/IF-1:200-1:500 |

### Target information (Q86UE4)


|                          |      |
|--------------------------|------|
| <b>Gene Name</b>         | MTDH |
| <b>About Gene Symbol</b> |      |




Immunocytochemistry analysis of LYRIC (green) in U87-MG using LYRIC antibody, and DAPI (blue).



Western blot analysis of LYRIC/AEG1 in 3T3, PC-12 lysates using LYRIC/AEG1 antibody.



Western blot analysis of LYRIC/AEG1 in CHO-K1, Ramos lysates using LYRIC/AEG1 antibody



Immunohistochemistry analysis of paraffin-embedded Human tonsil using LYRIC/AEG1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.