

## B7H4 Rabbit mAb


ABZ2014

### Product Overview

|                             |                                                            |
|-----------------------------|------------------------------------------------------------|
| <b>Name</b>                 | B7H4 Rabbit mAb                                            |
| <b>Accession(Primary)</b>   | Q7Z7D3                                                     |
| <b>Host Species</b>         | Rabbit                                                     |
| <b>Clonality</b>            | Monoclonal                                                 |
| <b>Isotype</b>              | IgG                                                        |
| <b>Species Reactivity</b>   | Human                                                      |
| <b>Calculated MW</b>        | Calculated MW: 31 kDa; Observed MW: 70 kDa                 |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:200-1:1000;<br>ICC/IF-1:100-1:200 |

### Target information (Q7Z7D3)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | VTCN1 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of B7H4 in MDA-MB-468,  
T47D lysates using B7H4 antibody.