

Aprataxin Rabbit mAb


ABZ2015

Product Overview


Name	Aprataxin Rabbit mAb
Accession(Primary)	Q7Z2E3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; ICC/IF-1:100-1:200

Target information (Q7Z2E3)

Gene Name	APTX
About Gene Symbol	



Immunocytochemistry analysis of Aprataxin (green) in U87-MG using Aprataxin antibody, and DAPI(blue).



Western blot analysis of Aprataxin in K562, HeLa lysates using Aprataxin antibody.