

CD272/BTLA Rabbit mAb


ABZ2016

Product Overview


Name	CD272/BTLA Rabbit mAb
Accession(Primary)	Q7TSA3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 34 kDa; Observed MW: 55-80 kDa
Recommended Dilution	WB-1:500-1:1000; ICC/IF-1:100; FC-1:200-1:500; IP-1:20-1:50; IHC-1:50-1:100

Target information (Q7TSA3)


Gene Name	Btla
About Gene Symbol	



Immunohistochemistry analysis of paraffin-embedded mouse spleen using CD272 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of CD272/BTLA in 3T3 lysates using CD272/BTLA antibody.



Western blot analysis of CD272/BTLA in K562,
HeLa lysates using CD272/BTLA antibody