


TAOK1 Rabbit mAb

ABZ2017


Product Overview

Name	TAOK1 Rabbit mAb
Accession(Primary)	Q7L7X3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 116 kDa; Observed MW: 116 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (Q7L7X3)

Gene Name TAOK1**About Gene Symbol**

Western blot analysis of TAOK1 in 3T3, Raw264.7 lysates using TAOK1 antibody.



Western blot analysis of TAOK1 in C6, Ramos
lysates using TAOK1 antibody.