

## CYFIP1 Rabbit mAb


ABZ2018

### Product Overview

|                             |                                                 |
|-----------------------------|-------------------------------------------------|
| <b>Name</b>                 | CYFIP1 Rabbit mAb                               |
| <b>Accession(Primary)</b>   | Q7L576                                          |
| <b>Host Species</b>         | Rabbit                                          |
| <b>Clonality</b>            | Monoclonal                                      |
| <b>Isotype</b>              | IgG                                             |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                               |
| <b>Calculated MW</b>        | Calculated MW: 145 kDa; Observed MW: 145 kDa    |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; ICC/IF-1:20-1:50; FC-1:20-1:50 |

### Target information (Q7L576)

|                          |                                                                                                                                                                                                                                                                                                                       |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | CYFIP1                                                                                                                                                                                                                                                                                                                |
| <b>About Gene Symbol</b> | CYFIP1 (Cytoplasmic FMR1 Interacting Protein 1) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of CYFIP1 in K562, rat Brain, C6, 3T3, HeLa lysates using CYFIP1 antibody.