

KCTD9 Rabbit mAb

ABZ2019


Product Overview

Name	KCTD9 Rabbit mAb
Accession(Primary)	Q7L273
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 43 kDa; Observed MW: 43 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:20-1:50


Target information (Q7L273)

Gene Name KCTD9

About Gene Symbol



Western blot analysis of KCTD9 in K562, HeLa lysates using KCTD9 antibody.



Immunocytochemistry analysis of KCTD9(green)
in HeLa using KCTD9 antibody, and DAPI(blue)