

RBM43 Rabbit mAb

ABZ2020


Product Overview

Name	RBM43 Rabbit mAb
Accession(Primary)	Q6ZSC3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q6ZSC3)

Gene Name RBM43

About Gene Symbol



Western blot analysis of RBM43 in HeLa lysates using RBM43 antibody.

