

# Nesprin3 Rabbit mAb


ABZ2021

## Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | Nesprin3 Rabbit mAb                          |
| <b>Accession(Primary)</b>   | Q6ZMZ3                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human                                        |
| <b>Calculated MW</b>        | Calculated MW: 112 kDa; Observed MW: 112 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                             |

## Target information (Q6ZMZ3)

|                          |                                                                                                                                                                                                                                                                                                                                           |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | SYNE3                                                                                                                                                                                                                                                                                                                                     |
| <b>About Gene Symbol</b> | SYNE3 (Spectrin Repeat Containing Nuclear Envelope Family Member 3) is a structural cytoskeletal protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Nesprin3 in Jurkat, HeLa lysates using Nesprin3 antibody.