

CRTC3 Rabbit mAb


ABZ2022

Product Overview


Name	CRTC3 Rabbit mAb
Accession(Primary)	Q6UUV7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 67 kDa; Observed MW: 78 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100

Target information (Q6UUV7)

Gene Name	CRTC3
About Gene Symbol	CRTC3 (CREB Regulated Transcription Coactivator 3) is a transcriptional coactivator involved in gene transcription and regulation. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate c



Immunocytochemistry analysis of CRTC3 (green) in HeLa using CRTC3 antibody, and DAPI (blue).



Western blot analysis of CRTC3 in K562, rat Brain, C6, 3T3, HeLa lysates using CRTC3 antibody.