

ULK3 Rabbit mAb


ABZ2024

Product Overview

Name	ULK3 Rabbit mAb
Accession(Primary)	Q6PHR2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:100-1:200; FC-1:10-1:100

Target information (Q6PHR2)

Gene Name	ULK3
About Gene Symbol	



Western blot analysis of ULK3 in K562, HeLa
lysates using ULK3 antibody.