

SLC27A4 / FATP4 Rabbit mAb


ABZ2025

Product Overview


Name	SLC27A4 / FATP4 Rabbit mAb
Accession(Primary)	Q6P1M0
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 72 kDa; Observed MW: 72 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:20-1:50; ICC/IF-1:50-1:100; IP-1:20-1:50

Target information (Q6P1M0)

Gene Name	SLC27A4
About Gene Symbol	



Western blot analysis of SLC27A4/FATP4 in K562, HeLa lysates using SLC27A4 antibody.



Immunocytochemistry analysis of SLC27A4/FATP4(green) in HeLa using SLC27A4/FATP4 antibody, and DAPI(blue)