

PMF-1 Rabbit mAb


ABZ2026

Product Overview


Name	PMF-1 Rabbit mAb
Accession(Primary)	Q6P1K2
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 23 kDa; Observed MW: 23 kDa
Recommended Dilution	WB-1:500-1:1000; IP-1:50-1:100

Target information (Q6P1K2)

Gene Name	PMF1
About Gene Symbol	PMF1 (Polyamine Modulated Factor 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PMF1 in C6, 3T3, Raw264.7 lysates using PMF 1 antibody.



Western blot analysis of PMF1 in HeLa, A549 lysates using PMF1 antibody.