

ALKBH2 Rabbit mAb

ABZ2027


Product Overview

Name	ALKBH2 Rabbit mAb
Accession(Primary)	Q6NS38
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 29 kDa; Observed MW: 29 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q6NS38)

Gene Name ALKBH2

About Gene Symbol



Western blot analysis of ALKBH2 in C6, 3T3, HeLa lysates using ALKBH2 antibody.

