

NCAPH2 Rabbit mAb


ABZ2028

Product Overview

Name	NCAPH2 Rabbit mAb
Accession(Primary)	Q6IBW4
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 68 kDa; Observed MW: 68 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200; IP-1:20

Target information (Q6IBW4)

Gene Name	NCAPH2
About Gene Symbol	



Western blot analysis of NCAPH2 in HeLa, Jurkat
lysates using NCAPH2 antibody.