


## N4BP1 Rabbit mAb

ABZ2029

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | N4BP1 Rabbit mAb                           |
| <b>Accession(Primary)</b>   | Q6A037                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Mouse, Rat                                 |
| <b>Calculated MW</b>        | Calculated MW: 99 kDa; Observed MW: 99 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; FC-1:100-1:500           |

### Target information (Q6A037)

**Gene Name** N4Bp1**About Gene Symbol**

Western blot analysis of N4BP1 in K562, HeLa lysates using N4BP1 antibody.

