


FADD Rabbit mAb

ABZ2031


Product Overview

Name	FADD Rabbit mAb
Accession(Primary)	Q61160
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse
Calculated MW	Calculated MW: 23 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (Q61160)

Gene Name Fadd**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded mouse liver using FADD antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of FADD in Raw264.7
lysates using FADD antibody.



Western blot analysis of FADD in SH-SY5Y lysates
using FADD antibody