

Adiponectin Rabbit mAb


ABZ2032

Product Overview


Name	Adiponectin Rabbit mAb
Accession(Primary)	Q60994
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 27 kDa; Observed MW: 30 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:500-1:2000; ICC/IF-1:500-1:1000; FC-1:20-1:50

Target information (Q60994)

Gene Name	Adipoq
About Gene Symbol	



Western blot analysis of Adiponectin in rat fat lysates using Adiponectin antibody.



Western blot analysis of Adiponectin in 3T3-L1 lysates using Adiponectin antibody.