

## XIAP Rabbit mAb

ABZ2033


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | XIAP Rabbit mAb                            |
| <b>Accession(Primary)</b>   | Q60989                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 56 kDa; Observed MW: 53 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

### Target information (Q60989)

**Gene Name** Xiap

**About Gene Symbol**



Western blot analysis of XIAP in K562, rat Brain, 3T3 lysates using XIAP antibody.

